


Programming within BrainVISA project






What is the BrainVISA project ?




- Complete development environment
- RedMine based forge
- Integration of various programming languages
- Multiplatform programming
- Documentation management
- Packaging and distribution

Typical structure of a BrainVISA project






- Part I – Introduction
- Part II – Programming with BrainVISA
- Part III – Programming with Anatomist
- Part IV – Programming with Aims in Python
- Part V – A complete example




Part I - Introduction







- User environnements
- How to start with Python ?
- Useful modules
- Exercises




- Shell
 - Mainly used to launch command lines
- Python
 - To be used only if IPython is not available
- **IPython**
 - Should be the main developement environnement
- BrainVISA
 - IPython plus BrainVISA specific extensions
- Text editor
 - Linux: kate or gedit



- Tutorial
 - <http://docs.python.org/tutorial/>
- Introduction to Python for Science
 - <http://gael-varoquaux.info/python4science.pdf>
- Some videos
 - <http://www.archive.org/search.php?query=SciPy%202009%20tutorial>
- In french
 - <http://dakarlug.org/pat/scientifique/html/index.html>



- Standard Python modules (not included by default)
 - SciPy, NumPy
 - Matplotlib
 - PyQt
- BrainVISA project modules
 - soma, aims
 - anatomist
 - brainvisa




- PBV_1: Count the number of occurrences of "BrainVISA" in web page <http://brainvisa.info>

```
import urllib
occurences = 0
file = urllib.urlopen( 'http://brainvisa.info/' )
for line in file:
    occurences = occurences + line.count( 'brainvisa' )
print occurences
```

- PBV_2: Show a diagram representing the size of BrainVISA bibliography pages between 2002 and 2009


```
import urllib
from pylab import *
url = 'http://brainvisa.info/biblio/en/Year/200%d.html'
years = range( 2, 10 )
values = [ len( urllib.urlopen( url % year ).read() ) for year in years ]
plot( years, values, linewidth=1.0)
ylabel('Size')
xlabel('Year')
title('Size of BrainVISA publication page per year')
grid(True)
show()
```






Part II – Programming with BrainVISA







- II.1 – Processes
- II.2 – Databases and Ontologies
- II.3 – Toolboxes




II.1 – Processes




BrainVISA main concepts



What is a process ?




- A process is an interface between BrainVISA and another software





- Reusable component. For instance, a process can call other processes.
- Python file containing one main function


How to create a new process ?




Create a file with **.py** extension in a toolbox **processes** directory




Special toolbox: My processes



- Located in user directory
 - Linux & Mac: `$HOME/.brainvisa`
 - Windows XP: `C:\Documents and Settings\<login>\.brainvisa`
 - Windows VISTA: `C:\Users\<login>\.brainvisa`
- Behave like a regular toolbox



Structure of a BrainVISA process



```
# -*- coding: utf-8 -*-  
from neuroProcesses import *
```

header

```
signature = Signature(  
    'image_input', ReadDiskItem(  
        '4D Volume',  
        [ 'NIFTI-1 image', 'SPM image', 'DICOM image', 'GIS image' ] ),  
    'image_output', WriteDiskItem(  
        '4D Volume',  
        'Aims writable volume formats' ),  
    'mode', Choice ( ( 'less than', 'lt' ),  
        ( 'less or equal', 'le' ),  
        ( 'greater than', 'gt' ),  
        ( 'greater or equal', 'ge' ),  
        ( 'equal', 'eq' ),  
        ( 'different', 'di' ),  
        ( 'between', 'be' ),  
        ( 'outside', 'ou' ) ),  
    'threshold1', Float(),  
    'threshold2', Float(),  
    'binary', Boolean(),  
)
```


signature


```
def initialization( self ):  
    self.setOptional( 'threshold2', 'binary' )  
    self.binary = 0  
    self.threshold1=0  
    self.mode='gt'
```

initialization

```
def execution( self, context ):  
    command = [ 'AimsThreshold',  
        '-i', self.image_input,  
        '-o', self.image_output,  
        '-m', self.mode,  
        '-t', self.threshold1 ]  
    if self.threshold2:  
        command += [ '-u', self.threshold2]  
    if self.binary:  
        command += [ '-b']  
    context.system( *command )
```

body






Example of a signature for a thresholding process

```
# -*- coding: utf-8 -*-  
from neuroProcesses import *  
name = 'Human readable error'  
userLevel = 0  
roles = ( 'viewer', )  
def validation( self ):  
    try:  
        import my_package.my_module  
    except:  
        raise ValidationError('my_module is not installed')
```


0 = Visible to everybody
1 = Invisible to basic users
2 = Invisible to basic and
advanced users

Possible roles:
viewer
editor
converter
importer

Signature of a BrainVISA process (1/2)




- Defines the parameters of a process
- `signature = Signature(parameter_list)`
 - `parameter_list` \leftarrow name, type, name, type, ...
 - Parameter types:
 - `String()`
 - `Number()`, `Float()`, `Integer()`
 - `Boolean()`
 - `Choice(<value>, (<label>, <value>), ...)`
 - `ReadDiskItem(type, formats)`
 - `WriteDiskItem(type, formats)`



Example of signature for a thresholding process

```
signature = Signature(  
    'image_input', ReadDiskItem(  
        '4D Volume',  
        [ 'NIFTI-1 image', 'SPM image', 'DICOM image', 'GIS image' ] ),  
    'image_output', WriteDiskItem(  
        '4D Volume',  
        'Aims writable volume formats' ),  
    'mode', Choice ( ( 'less than', 'lt' ),  
        ( 'less or equal', 'le' ),  
        ( 'greater than', 'gt' ),  
        ( 'greater or equal', 'ge' ),  
        ( 'equal', 'eq' ),  
        ( 'different', 'di' ),  
        ( 'between', 'be' ),  
        ( 'outside', 'ou' ) ),  
    'threshold1', Float(),  
    'threshold2', Float(),  
    'binary', Boolean(),  
)
```



Initialization example

```
def initialization( self ):  
    self.binary = 0  
    self.threshold1=0  
    self.mode='gt'  
    self.setOptional( 'threshold2', 'binary' )  
    self.linkParameters( 'threshold2', 'threshold1' )  
  
    eNode = SerialExecutionNode( self.name,  
        parameterized=self )  
    eNode.addChild( 'correction',  
        ProcessExecutionNode( 'DiffusionEPICorrection',  
                                optional=True, selected=True ) )  
    eNode.addChild( 'create_mask',  
        ProcessExecutionNode( 'DiffusionT2BrainMask',  
                                optional=True, selected=True ) )
```


Parameters initialization

Optional parameters

Linked parameters

Pipeline
creation

The real work : the process body



Python function


```
def execution( self, context ):
    command = [ 'AimsThreshold',
                '-i', self.image_input,
                '-o', self.image_output,
                '-m', self.mode,
                '-t', self.threshold1 ]
    if self.threshold2:
        command += [ '-u', self.threshold2 ]
    if self.binary:
        command += [ '-b' ]
    context.system( *command )
```

self : get parameters values


Execution context:

Manage interactions with system and user in a controlled and context dependent way


- write, warning, error : prints a message.
- log : writes a message in the BrainVISA log file.
- ask, dialog : asks a question to the user.
- temporary : creates a temporary file.
- system: call a system command.
- runProcess : runs a process.
- matlab : calls a Matlab command.




- PBV_3: create a process taking an **image** and a **threshold** and writing a **binary image** containing all voxels **greater than** threshold.
- PBV_4: create a process that takes a **binary image** and create a **3D mesh** from this image.
- PBV_5: create a process **combining PBV_3 and PBV_4**. Input is image and threshold, output is mesh. Thresholded image is only used internally in a temporary file.
- PBV_6: create a process that **display** an **image** and a **mesh** in the same Anatomist window



- Pipeline: combination of existing processes
- Customized parameter links
- Dynamic signature
- Customization of graphical interface




- PBV_7: create a **pipeline** process **chaining PBV_3 and PBV_4**.



II.2 – Databases and Ontologies




Why using BrainVISA databases ?




- Define a shared data organization
- Reuse data from people who left the lab
- Share data with people from other labs
- Make links between data
 - Find the head mesh corresponding to this MRI scan
- Automation of data processing


What is a BrainVISA database ?




- What is a BrainVISA database ?
 - A **directory** containing **data files** that are organized in a **hierarchy** that follows an **ontology**
 - A **relational database** built from an **ontology** and allowing to make efficient **selection** requests on **data files**.
- Two ways of seeing the same data connected together by a common ontology




bijection



What is a BrainVISA ontology ?



Hierarchy of types



Attributes

Raw T1 MRI
protocol
subject
acquisition
normalization

T1 MRI Bias Corrected:
protocol
subject
acquisition
analysis

T1 MRI:
protocol
subject
acquisition
normalization
analysis

Data selection from databases rely only on ontology

Select data from databases



Resultats du filtre / Nom(s) suggere(s)
6 item(s)

	type	protocol	subject	acquisition	normalization	analysis	format	database
1	Raw T1 MRI	demo	sujet01	200810			GIS image	/home/yann/jirfni2008/db_processed
2	Raw T1 MRI	demo	sujet02	200810			GIS image	/home/yann/jirfni2008/db_processed
3	Raw T1 MRI	demo	sujet03	200810			GIS image	/home/yann/jirfni2008/db_processed
4	T1 MRI Bias Corrected	demo	sujet01	200810		default_analysis	GIS image	/home/yann/jirfni2008/db_processed
5	T1 MRI Bias Corrected	demo	sujet02	200810		default_analysis	GIS image	/home/yann/jirfni2008/db_processed
6	T1 MRI Bias Corrected	demo	sujet03	200810		default_analysis	GIS image	/home/yann/jirfni2008/db_processed

infos

Ok Cancel

BrainVISA database hierarchy




Raw T1 MRI:

{protocol}/{subject}/t1mri/{acquisition}/<subject>

{protocol}/{subject}/t1mri/{acquisition}/normalized_{normalization}_<subject>

T1 MRI Bias Corrected:

{protocol}/{subject}/t1mri/{acquisition}/{analysis}/nobias_<subject>




Create types in <toolbox>/types/*.py


```
include( 'builtin' )  
include( 'anatomy' )  
  
FileType( 'T1 MRI Bias Corrected', 'T1 MRI' )  
  
Format( 'BrainVISA Gyri Model', "f|*.gyr" )  
FileType( 'Gyri Model', 'Any Type', BrainVISA Gyri Model' )
```

Add hierarchy rules in <toolbox>/hierarchies/brainvisa-3.1.0/*.py

```
include( 'base' )  
  
insert( {protocol}/{subject}/t1mri/{acquisition}/{analysis}',  
        'nobias_<subject>', SetType( 'T1 MRI Bias Corrected' ),  
        )
```


- PBV_8: Create a **new type**: "Mesh from Threshold"
- PBV_9: Create a **new hierarchy entry** for "Mesh from Threshold" with the same key attributes as "Fractional Anisotropy"
- PBV_10: Create a **process** that create a mesh of type "**Mesh from Threshold**" from a thresholded "Fractional Anisotropy" image. Include a **link** between input image and output mesh.
- PBV_11: Create a **viewer** for "Mesh from threshold" that display the mesh with the corresponding FA image.



II.3 – Toolboxes




What is a BrainVISA toolbox ?




- An set of BrainVISA extensions
- Processes
- Ontology
- Documentation


Why creating a toolbox ?




- Organize processes according to a topic
- Distribute BrainVISA extensions




- BrainVISA is in CeCill v2 (i.e. GPL)
- Therefore a BrainVISA toolbox must be in GPL
- However, processing libraries can have any licence as long as they do not rely on BrainVISA and can be distributed separately
- Example: T1 MRI Segementation Toolbox







- Create a new **<name>** directory in **toolboxes** directory
- Put processes in **processes** directory
- Put ontology extension in **types** and **hierarchies** directories
- Create a configuration file: **<name>.py**
 - `userName = 'User will see this name'`
 - `icon = file path of the icon that will represent the toolbox in graphical interface. Optional, there is a default icon.`
 - `description` : tooltip for the toolbox (default is the name of the toolbox)




- The same process can appear in several toolboxes
- Links are stored in a minf file: **<name>.minf**
- Example :
 - Diffusion importation processes from Diffusion & tracking toolbox
 - These processes are also in Data management toolbox
- BrainVISA interface can be used to create this file



- Structure of BrainVISA package directory:
 - bin : executables
 - brainvisa/toolboxes : BrainVISA toolboxes
 - include : C/C++ headers
 - lib : dynamic libraries
 - python : Python libraries
 - share : shared data and documentation
- Create an archive that adds files to the existing structure




- PBV_12: Create a package containing BrainVISA extensions corresponding to exercises PBV_8, PBV_9, PBV_10 and PBV_11



Part III – Programming with Anatomist



Anatomist tutorial – programming part




Tutorial about Anatomist python API:

http://brainvisa.info/doc/anatomist/ana_training/en/html/ch08.html




Part IV – Programming with Aims in Python






Tutorial about AIMS python API:

http://brainvisa.info/doc/aimsdata/aims_training/en/html/ch10.html



Part V – A complete example





- Part 1: Create a BrainVISA process to compute the coordinates of the maximum using aims


```
from neuroProcesses import *
import neuroConfig
from brainvisa import anatomist
import numpy

signature = Signature(
    'image', ReadDiskItem( '4D Volume', 'Aims readable volume formats' ),
)

def execution( self, context ):

    # Read input image
    aims_image = aims.read( self.image.fullPath() )

    # Compute coordinate of maximum voxel with Python
    aims_max_coordinate = ( 0, 0, 0, 0 )
    maximum = image.at( 0, 0, 0, 0 )
    for t in xrange( image.getSizeT() ):
        for z in xrange( image.getSizeZ() ):
            for y in xrange( image.getSizeY() ):
                for x in xrange( image.getSizeX() ):
                    v = aims_image.at( x, y, z, t )
                    if v > maximum:
                        maximum = v
                        aims_max_coordinate = ( x, y, z, t )
    context.write( 'Aims says max is', maximum, 'at', aims_max_coordinate )
```



- Part 2: Use numpy to process maximum coordinates


```
from neuroProcesses import *
import neuroConfig
from brainvisa import anatomist
import numpy

signature = Signature(
    'image', ReadDiskItem( '4D Volume', 'Aims readable volume formats' ),
)

def execution( self, context ):

    # Read input image
    aims_image = aims.read( self.image.fullPath() )

    numpy_matrix = numpy.array(aims_image, copy = False )
    numpy_max_coordinate = numpy.unravel_index( numpy_matrix.argmax(), numpy_matrix.shape )
    context.write( 'Numpy says max is', aims_image.at( *numpy_max_coordinate ), 'at', numpy_max_coordinate )
```


- Part 3: Add read of a spherical mesh of size 100 located in standard BrainVISA directory and write it in a temporary file

```
aims_sphere = aims.read( os.path.join( neuroConfig.dataPath[0].directory, 'standardmeshes', 'ico100_7.mesh') )

# Change the size of the sphere and center it on image maximum
voxel_size = aims_image.header()[ 'voxel_size' ]
for vertex in aims_sphere.vertex():
    vertex *= self.sphere_size / 100.0
    vertex[ 0 ] += aims_max_coordinate[ 0 ] * voxel_size[ 0 ]
    vertex[ 1 ] += aims_max_coordinate[ 1 ] * voxel_size[ 1 ]
    vertex[ 2 ] += aims_max_coordinate[ 2 ] * voxel_size[ 2 ]

# Create a temporary file name for the modified sphere
diskitem_sphere = context.temporary( 'Mesh Mesh' )

# Write the modified sphere in temporary file
aims.write( aims_sphere, diskitem_sphere.fullPath() )
```




- Part 4: Create a texture on the sphere that takes the value from an interpolator for each point of the sphere and save it in a temporary file

```
# Create an interpolator on the image to be able to get a value for any millimeter coordinate
interpolator = aims.getLinearInterpolator( aims_image )


# Create a texture on the sphere that takes the value from the interpolator for each point of the sphere
vertices_count = aims_sphere.vertex().size()
aims_texture = aims.Texture_FLOAT( vertices_count )
for i in xrange( vertices_count ):
    aims_texture[ i ] = interpolator.value( *aims_sphere.vertex()[ i ] )

# Write the created texture to a temporary file
diskitem_texture = context.temporary( 'Texture' )
aims_time_texture = aims.TimeTexture_FLOAT()
aims_time_texture[ 0 ] = aims_texture
aims.write( aims_time_texture, diskitem_texture.fullPath() )
```



- Part 5: Open anatomist and visualize the textured sphere at the maximum of the image

```
# Open Anatomist
a = anatomist.Anatomist()
# Load sphere in anatomist from temporary file
aSphere = a.loadObject( diskitem_sphere.fullPath() )
# Load sphere texture in anatomist from temporary file
aTexture = a.loadObject( diskitem_texture.fullPath() )
# Fusion sphere and texture to create a textured object
aTexturedSphere = a.fusionObjects( (aSphere, aTexture), 'FusionTexSurfMethod' )
# Load image in anatomist
aImage = a.loadObject( self.image.fullPath() )
# Create an Axial window
aWindow = a.createWindow( 'Axial' )
# Display textured sphere and image in window
aWindow.addObject( (aTexturedSphere, aImage) )
# Move Anatomist cursor to the center of the voxel with maximum value
aWindow.moveLinkedCursor( ( aims_max_coordinate[ 0 ] * voxel_size[ 0 ],
                           aims_max_coordinate[ 1 ] * voxel_size[ 1 ],
                           aims_max_coordinate[ 2 ] * voxel_size[ 2 ] ) )
```



- Part 6: finalize the process to not destroy python objects used

```
# Return objects that must not be destroyed immediately  
return [ aTexturedSphere, almage, aWindow ]
```